NEDERLANDS

Beppy Pure® tampon

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ENGLISH

Beppy Pure® tampon

Comfortable tampon without string.

The Beppy Pure® tampon is developed for women who want to feel free during their menstruation. The Beppy Pure® tampon is a cordless sponge tampon. Therefore the tampon is ideal for use during a sauna visit, sports, swimming, sexual intercourse and more. The shape of the tampon adjusts to your body due to its soft and flexible structure.

STEP 1: Wash your hands always before and after inserting the tampon. Open the foil at the notch.

STEP 2: Fold the tampon between your thumb and index finger. Make sure that the integrated grip in the tampon is pointed downwards.

STEP 3 – inserting the tampon: Adopt a comfortable position; the tampon can be inserted while standing with one leg raised above the other, crouching or while sitting. Relax the pelvic floor muscles and insert the tampon. Push the tampon with your middle or index finger in the direction of your tailbone until it reaches the cervix.

STEP 4 – removing the tampon: Wash your hands always before and after removing the tampon. Relax and take a crouching position. Use your index and middle finger like tweezers to grab the tampon and remove the tampon slowly. You can also use the integrated grip in the tampon. Press slightly when removing the tampon.

How long can you use the tampon?

It is different per woman when the tampon should be removed/replaced. It depends on the amount of blood loss and the day of your menstrual period. Guideline: every 2-4 hours during heavy flow, every 4-8 hours during the days of lighter flow.

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
BP10121- V1922

All rights reserved © Asha International. Designed and made in Holland (EU). Registered brand names and design. Last updated May 2019.

Gooi de tampon alsjeblieft in de prullenmand / Please, throw the tampon in the trash bin / Bitte werfen Sie den Tampon in den Mülleimer / Prière de jeter le tampon dans la poubelle / Por favor, arroja el tampón a la papelera / Após o uso, não jogue o absorvente Beppy no vaso sanitário.

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ENGLISH

Please note: If you leave the tampon in the vagina for more than 8 hours, its safety and effectiveness is no longer guaranteed. The tampon can be used only once. Reusing them may cause infection. The tampon is not a contraceptive. Do not use tampons between your periods for regular discharge. Do not use when your hymen is intact and/or when the use-by date on the pack has expired.

IMPORTANT about Toxic Shock Syndrome (TSS):

TSS is a rare but serious and sometimes fatal disease, which is associated with the use of menstrual tampons.

What are the symptoms of TSS? Early symptoms may begin suddenly and are similar to the flu: a sudden fever - usually 39°C or higher, vomiting, diarrhoea, fainting, dizziness or a rash that looks like sunburn. Early recognition and treatment is very important and can generally prevent serious illness.

What should you do if you think you have TSS? If, during your period or shortly thereafter, you have any of the above symptoms, **REMOVE YOUR TAMPON AND CONSULT A DOCTOR IMMEDIATELY.** It is very important to tell your doctor that you have your period and have used tampons.

Who is at risk? TSS can occur in all women who use tampons during their periods. For women under 30, higher risks have been reported.

Reducing or avoiding the risk? To reduce the risk of TSS, use tampons with the lowest absorbency to suit your flow. To avoid or reduce the risk, do not use tampons or alternate the use of tampons with sanitary towels.

Consult a doctor if you want to use tampons again and have had symptoms of TSS in the past or have questions about TSS. As a general rule, it is necessary to visit your doctor, gynaecologist or A&E if you cannot remove a tampon.

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